




**FEMALE PATIENT CARE**

# Biopsy Punches Colposcopy Sets





## WHY PHYSICIANS CHOOSE SKLAR INSTRUMENTS:

### Quality, Service, Price and Sklar Experience

Gynecologists know that choosing quality surgical instruments takes more than flipping through a catalog and placing an order. While price is important, there is no substitution for quality. Sklar gynecological instruments are overwhelmingly chosen by discerning surgeons because of Sklar's high quality standards, unparalleled service, more than 100 years of experience, and competitive pricing.

**QUALITY** - All Sklar instruments are guaranteed for life against manufacturing defects. Our instruments are carefully examined in the factory using a 51-point checklist: if an instrument does not meet Sklar standards (which exceed industry standards), it does not leave the factory.

**SERVICE** - Our customer service is unequalled in the surgical instrument industry. Our thoroughly-trained and knowledgeable team of professionals is ready to provide each customer with personalized attention to answer questions and meet your needs.

**EXPERIENCE** - With more than 100 years experience in the research, design, and manufacture of quality surgical instruments, Sklar offers advanced technical know-how in all facets of surgical instrumentation.

**PRICE** - That same 100 years of experience has enabled Sklar to produce superior products that are affordable.

With Sklar precision instruments, you concentrate on your patient and surgery, not the instrument you're using.





## FEMALE PATIENT CARE

### Biopsy Punches and Colposcopy Sets

**Precision engineered,** hand-crafted and imported from Germany, Sklar Biopsy Punches and Colposcopy Instruments enjoy a reputation as being the finest in the field. A complete selection of forcep tips combined with three specific handle configurations provide the flexibility, comfort and control required by the most selective physicians.

In addition to three handle configurations, Sklar also offers a line of fully rotational biopsy forceps designed to obtain a perfect sample while not obstructing the biopsy site. The shaft is detachable from the handle, rotates a full 360°, and is offered in a full assortment of tip-shaft combinations.

**SKLAR OR GRADE** biopsy punches insure the highest level of hand-crafted workmanship. The precision, balance and feel of these fine instruments is apparent from the moment the physician lifts it from the tray.

**MICRO-CUT PREMIUM GRADE** biopsy punches combine the same fine craftsmanship of our OR Grade with a knife-honed edge on the upper jaw for a cleaner and more precise cut.

**MICRO-CUT**

As with all our fine instruments, Sklar combines over a century of experience with precision craftsmanship and highly competitive pricing to provide the finest value in quality gynecological instruments.

## SKLAR BIOPSY FORCEPS

	<b>TISCHLER-KEVORKIAN</b>	3mm x 9.5mm Bite	■
	<b>TISCHLER-KEVORKIAN</b>	3mm x 9.5mm Wide Bite	■
	<b>MINI TISCHLER</b>	2.3mm x 4.2mm Bite	●
	<b>MINI TISCHLER</b>	2.3mm x 4.2mm Wide Bite	●
	<b>TISCHLER OVAL</b>	3mm x 7mm Bite	●
	<b>TISCHLER OVAL</b>	3mm x 7mm Wide Bite	●
	<b>KEVORKIAN</b>	3mm x 8mm Bite	■
	<b>KEVORKIAN</b>	3mm x 8mm Wide Bite	■

### SKLAR BIOPSY FORCEPS

#### OR GRADE

- 90-6887** Tischler-Kevorkian Biopsy Forcep
- 90-9570** Tischler-Kevorkian Biopsy Forcep, Wide Bite
- 90-7530** Mini Tischler Biopsy Forcep
- 90-9555** Mini Tischler Biopsy Forcep, Wide Bite
- 90-7520** Tischler Oval Biopsy Forcep
- 90-9560** Tischler Oval Biopsy Forcep, Wide Bite
- 90-7540** Kevorkian Biopsy Forcep
- 90-9575** Kevorkian Biopsy Forcep, Wide Bite

#### PREMIUM GRADE

#### MICRO-CUT

- 90-8920** Micro-Cut Tischler-Kevorkian Biopsy Forcep
- 90-8945** Micro-Cut Mini Tischler Biopsy Forcep
- 90-8940** Micro-Cut Tischler Oval Biopsy Forcep
- 90-8952** Micro-Cut Kevorkian Biopsy Forcep



# SKLAR BIOPSY FORCEPS



## OVAL

3mm x 5mm Bite



## OVAL

3mm x 5mm Wide Bite



## OVAL BIOPSY FORCEPS

### OR GRADE

**90-7525** Oval Biopsy Forcep

**90-9565** Oval Biopsy Forcep, Wide Bite

### PREMIUM GRADE

### MICRO-CUT

**90-8948** Micro-Cut Oval Biopsy Forcep

# SKLAR BIOPSY FORCEPS

	<b>EPPENDORFER</b>	
	4mm x 8mm Bite	
	<b>KEVORKIAN-YOUNGE</b>	
	3mm x 9.7mm Bite	
	<b>WITTNER STRAIGHT</b>	
	4mm x 8mm Bite	
	<b>WITTNER ANGULAR</b>	
	4mm x 8mm Bite	

## SKLAR BIOPSY FORCEPS

### OR GRADE

- 90-6820** Eppendorfer Biopsy Forcep
- 90-6610** Kevorkian-Younger Biopsy Forcep
- 90-6600** Wittner Biopsy Forcep - Straight
- 90-6602** Wittner Biopsy Forcep - Angular


### PREMIUM GRADE

### MICRO-CUT

- 90-8010** Micro-Cut Eppendorfer Biopsy Forcep
- 90-8905** Micro-Cut Kevorkian-Younger Biopsy Forcep
- 90-8900** Micro-Cut Wittner Biopsy Forcep - Straight



# SKLAR FULLY ROTATIONAL BIOPSY FORCEPS



- n 360° ROTATING SHAFT
- n COMPLETELY UNOBSTRUCTED BIOPSY SITE
- n SIMPLE, TIME-SAVING DETACHABLE SHAFT
- n FULL ASSORTMENT OF TIP-SHAFT CHOICES


<b>TISCHLER OVAL</b>		
3mm x 7mm Bite		
<b>OVAL</b>		
3mm x 5mm Bite		
<b>MINI TISCHLER</b>		
2.3mm x 4.2mm Bite		
<b>KEVORKIAN</b>		
3mm x 8mm Bite		

## FULLY ROTATIONAL BIOPSY FORCEPS

OR GRADE	
<b>90-7571</b>	Tischler Oval Biopsy Forcep w/ Handle
<b>90-7572</b>	Oval Biopsy Forcep w/ Handle
<b>90-7573</b>	Mini Tischler Biopsy Forcep w/ Handle
<b>90-7574</b>	Kevorkian Biopsy Forcep w/ Handle
<b>Cannula with Tip Only</b>	
<b>90-7576</b>	Tischler Oval
<b>90-7577</b>	Oval
<b>90-7578</b>	Mini Tischler
<b>90-7579</b>	Kevorkian
<b>90-7575</b>	Universal Handle Only

PREMIUM GRADE		MICRO-CUT
<b>90-8981</b>	Tischler Oval Biopsy Forcep w/ Handle	
<b>90-8982</b>	Oval Biopsy Forcep w/ Handle	
<b>90-8983</b>	Mini Tischler Biopsy Forcep w/ Handle	
<b>90-8984</b>	Kevorkian Biopsy Forcep w/ Handle	
<b>Cannula with Tip Only</b>		
<b>90-8986</b>	Tischler Oval	
<b>90-8987</b>	Oval	
<b>90-8988</b>	Mini Tischler	
<b>90-8989</b>	Kevorkian	
<b>90-8985</b>	Universal Gold Handle Only	

## COLPOSCOPY SETS



### STANDARD COLPOSCOPY SET

**Includes the following:**

**90-8557** Standard Colposcopy Set

**90-3995** Kogan Endocervical Speculum - Wide Wire Blades


**90-8334** Emmett (Iris) Hook - 4mm Diameter, 9" long, quick release

**90-6612** Kevorkian-Younger Endocervical Curette - 12" with Basket

**90-7520** Tischler Oval Biopsy Forcep - 3mm x 7mm Bite



## COLPOSCOPY SETS



### MINI COLPOSCOPY SET

#### Includes the following:

- 90-8555** Mini Colposcopy Set
- 90-3990** Kogan Endocervical Speculum - Narrow Wire Blades
- 90-8336** Iris Hook - 1mm Diameter, 9" long, quick release
- 90-6285** Mini Endocervical Curette
- 90-7530** Mini Tischler Biopsy Forcep - 2.3mm x 4.2mm Bite

## INSTRUMENT CARE AND CLEANING



### SKLAR INSTRUMENT CARE AND CLEANING PRODUCTS

- 10-1603** Sklar Soak, Instrument Decontaminate - Case of 4, 1 Gallon Bottles
- 10-2777** Sklar Enzymatic, Instrument Soak/Detergent - Case of 4, 1 Gallon Bottles
- 10-2774** SklarZyme, Enzymatic Detergent Concentrate - 1 Box of 32, 1/2 oz. Packets
- 10-1612** Sklar Kleen Liquid Concentrate - Case of 6, 8 oz. Bottles
- 10-1613** Sklar Kleen Liquid Concentrate - Case of 4, 1 Gallon Bottles
- 10-1630** Sklar Kleen Powder - Case of 6, 3 1/4 lb. Pitchers
- 10-1607** Sklar Kleen Powder - 5 lb. Pail
- 10-2701** Sklar Kleen LF Liquid Concentrate - Case of 4, 1 Gallon Bottles
- 10-1636** Sklar Lube - Case of 12 - 8 oz. Ready-to-use Spray Bottles
- 10-1635** Sklar Lube Concentrate - Case of 4, 1 Gallon Bottles
- 10-1650** Instrument Cleaning Brush - Stainless Steel - Pack of 3
- 10-1652** Instrument Cleaning Brush - Stainless Steel - Pack of 12

- 10-1651** Instrument Cleaning Brush - Nylon - Pack of 3
- 10-1657** Instrument Cleaning Brush - Nylon - Pack of 12
- 10-1626** Sklar Polish - Case of 12, 8 oz. Jars
- 10-1629** Sklar Polish - Case of 6, 8 oz. Jars
- 10-2762** Sklar Disinfectant - Case of 12, 24 oz. Spray Bottles
- 10-1643** Sklar Disinfectant - Case of 6, 24 oz. Spray Bottles
- 10-1653** Sklar Disinfectant Liquid - Case of 4, 1 Gallon Bottles
- 10-2800** Sklar Instru-Guard, 1-Step Detergent/Lubricant - Case of 4, 1 Gallon Bottles
- 10-1744** Sklar Instrument Tray - small
- 10-1724** Sklar Instrument Tray Cover - small
- 10-1746** Sklar Instrument Tray - large
- 10-1882** Sklar Instrument Tray Cover - large

## INSTRUMENT CLEANING PROCEDURE

### A. Sklar Soak - Disinfectant and Decontaminant

The use of *Sklar Soak* for 10 minutes substantially reduces the risk of exposure to many infectious organisms including TB and viruses including HIV-1 (Aids Virus), bacteria and fungi.

- Non-caustic formula
- Effective against HIV-1 (Aids Virus), TB, bacteria, and fungi
- EPA registered
- Meets all OSHA Standards on blood borne pathogen destruction
- Non-Glutaraldehyde!

### B. Cleaning

Sklar recognizes that user requirements vary. Therefore, Sklar now offers four custom detergent formulations.

- Sklar Enzymatic
- Sklar Kleen Liquid
- Sklar Kleen Powder
- Sklar Kleen LF (Low Foam)

### C. Sklar Instru-Guard Lube - Extends Instrument Life and Performance

Sklar Instru-Guard Lube is an easy to use, water soluble, non-silicone lubricant and rust inhibitor. When used daily as a second rinse after every cleaning, it keeps box locks operating smoothly, is effective even after autoclaving, and leaves no oily residue.

*Sklar's complete line of instrument care and cleaning products provides everything you need to keep your instruments performing like new.*

## INSTRUMENT CARE AND CLEANING PRODUCTS

### Sklar Enzymatic

*Sklar Enzymatic Detergent* rapidly dissolves organic debris on instruments and clothes. It is a low-foaming, quick cleaning detergent specifically designed for use in areas where instruments and/or cloth materials may become soiled with organic debris.

### Sklar Kleen

Specifically formulated for both manual and ultrasonic cleaning, *Sklar Kleen* provides superior wetting power to penetrate locks, serrations, and crevices of surgical instruments, glassware, plastics and rubber.

- Sklar Kleen Powder
- Sklar Kleen Liquid

### Sklar Polish

*Sklar Polish* is safe for all your stainless steel, chrome, nickel, brass, instrument trays, basins and racks. It also works well on porcelain and tile surfaces. The true value in regular *Sklar Polish* use is that it removes minerals from your instruments and other metal surfaces without changing the surface qualities of the items.

### Sklar Disinfectant

- Specially formulated for the control of cross contamination
- Effective against TB, Hydrophilic and Lipophilic viruses including HIV-1 (Aids Virus), bacteria and fungi on inanimate objects
- Safe non-toxic and non-irritating formula
- EPA registered
- Glutaraldehyde Free!



889 South Matlack Street  
West Chester, PA 19382

ph: 800.221.2166 / 610.430.3200  
fx: 610.696.9007  
e-mail: [international@sklarcorp.com](mailto:international@sklarcorp.com)  
[surgi@sklarcorp.com](mailto:surgi@sklarcorp.com)

[www.sklarcorp.com](http://www.sklarcorp.com)

