


**BAND BAGS  
SEE PAGE 11 FOR  
COMPLETE LISTING**



**SKLAR EQUIPMENT COVERS AND BAGS**

**OR PROTECTION AND  
HOUSEKEEPING SOLUTIONS**



# SKLAR DISINFECTANT™

Ready-To-Use  
Spray Surface Disinfectant

**2 Minute Kill Time!**

- Ready to use, non-toxic formula
- Specially formulated for the control of cross contamination
- Effective against a broad range of pathogens including H1N1, HIV-1, HPV, VRE, Norovirus, Hep B & C as well as TB, MRSA, A. Baumannii, and C. Diff (veg.).
- Safe for all areas of the facility including the NICU, operating rooms, isolation rooms, patient care areas and laboratories
- Satisfies OSHA's Bloodborne Pathogen Standards
- Non-irritating formula
- EPA registered

## PRODUCT NUMBER

10-1643

FEATURES: 24 OZ. SPRAY BOTTLES

QUANTITY: 6 / CASE

## PRODUCT NUMBER

10-1642

FEATURES: 24 OZ. SPRAY BOTTLE

QUANTITY: SINGLE

## PRODUCT NUMBER

10-1653

FEATURES: 1 GALLON BOTTLES

QUANTITY: 4 / CASE

## PRODUCT NUMBER

FEATURES: 1 GALLON BOTTLE

QUANTITY: SINGLE

## PRODUCT NUMBER

10-1644

FEATURES: 2.5 GALLON BOTTLES

QUANTITY: 2 / CASE



READY-TO-USE SPRAY SURFACE DISINFECTANT

**DISINFECTANT™**



**STERILIZATION POUCHES ..... 2**

**BIOHAZARD WASTE BAGS..... 3**

**LINEN BAGS ..... 4**

**SPECIMEN TRANSPORT BAGS..... 5**

**RECLOSABLE BAGS ..... 5**

**POST MORTEM BAGS ..... 6**

**FLUOROSCOPE COVERS..... 7**

**SETUP COVERS ..... 8**

**DOME COVERS ..... 9**

**MAYO STAND COVERS ..... 10**

**BAND BAGS ..... 11**

**WASTE CAN LINERS ..... 12**



©COPYRIGHT 2010. SKLAR INSTRUMENTS, WEST CHESTER, PA.  
 ALL RIGHTS RESERVED. NO PORTION OF THIS CATALOG MAY BE  
 REPRODUCED WITHOUT THE WRITTEN PERMISSION OF SKLAR.

# STERILIZATION POUCHES

PRODUCT NUMBER	FEATURES	SIZE	QUANTITY
10-6230	Self-Seal, Paper / Poly, Clear / White	3.5" x 9"	2000 Flat / Cs
10-6235	Self-Seal, Paper / Poly, Clear / White	3.5" x 9"	200 Flat / Cs
10-6240	Self-Seal, Paper / Poly, Clear / White	3.5" x 22"	1000 Flat / Cs
10-6245	Self-Seal, Paper / Poly, Clear / White	3.5" x 22"	200 Flat / Cs
10-6250	Self-Seal, Paper / Poly, Clear / White	5.25" x 10"	2000 Flat / Cs
10-6255	Self-Seal, Paper / Poly, Clear / White	5.25" x 10"	200 Flat / Cs
10-6260	Self-Seal, Paper / Poly, Clear / White	5" x 15"	1000 Flat / Cs
10-6265	Self-Seal, Paper / Poly, Clear / White	5" x 15"	200 Flat / Cs
10-6270	Self-Seal, Paper / Poly, Clear / White	7.5" x 13"	1000 Flat / Cs
10-6275	Self-Seal, Paper / Poly, Clear / White	7.5" x 13"	200 Flat / Cs
10-6280	Self-Seal, Paper / Poly, Clear / White	12" x 15"	1000 Flat / Cs
10-6285	Self-Seal, Paper / Poly, Clear / White	12" x 15"	200 Flat / Cs
10-6290	Self-Seal, Paper / Poly, Clear / White	12" x 18"	1000 Flat / Cs
10-6295	Self-Seal, Paper / Poly, Clear / White	12" x 18"	200 Flat / Cs
10-6300	Heat-Seal, Paper / Poly, Clear / White	4" x 8"	1000 Flat / Cs
10-6305	Heat-Seal, Paper / Poly, Clear / White	4" x 8"	200 Flat / Cs
10-6310	Heat-Seal, Paper / Poly, Clear / White	6" x 10"	1000 Flat / Cs
10-6315	Heat-Seal, Paper / Poly, Clear / White	6" x 10"	200 Flat / Cs
10-6320	Heat-Seal, Paper / Poly, Clear / White	6" x 15"	1000 Flat / Cs
10-6325	Heat-Seal, Paper / Poly, Clear / White	6" x 15"	200 Flat / Cs
10-6330	Heat-Seal, Paper / Poly, Clear / White	8" x 13"	1000 Flat / Cs
10-6335	Heat-Seal, Paper / Poly, Clear / White	8" x 13"	200 Flat / Cs
10-6340	Heat-Seal, Paper / Poly, Clear / White	12" x 15"	1000 Flat / Cs
10-6345	Heat-Seal, Paper / Poly, Clear / White	12" x 15"	200 Flat / Cs
10-6350	Self-Seal, Paper / Poly, Clear / White	3" x 100'	8 Rolls / Cs
10-6355	Self-Seal, Paper / Poly, Clear / White	4" x 100'	16 Rolls / Cs
10-6360	Self-Seal, Paper / Poly, Clear / White	6" x 100'	8 Rolls / Cs
10-6365	Self-Seal, Paper / Poly, Clear / White	10" x 100'	4 Rolls / Cs

Available in Self-Seal and Heat-Seal designs, Sklar Sterilization Pouches are made of high-quality, Medical Grade paper that resists fiber tears, triple rail seals to protect against penetration and allow for uniform peel when opening, and features a see-through front to allow visual inspection of contents. Easy-to-read color changing process indicators are compatible with Steam or ETO sterilization.

## CROSS REFERENCE

SKLAR	MEDICAL ACTION	MEDLINE	HEALTHMARK	CARDINAL
10-6230	520	MPP100520	SS-T1	92309
10-6235	520	MPP100520Z	SS-T1	92309
10-6240	565	MPP100565	SS-T3	92322
10-6245	565	MPP100565Z	SS-T3	92322
10-6250	525	MPP100525	SS-T4	92510
10-6255	525	MPP100525Z	SS-T4	92510
10-6260	530	MPP100550	SS-T4A	92515
10-6265	530	MPP100550Z	SS-T4A	92515
10-6270	545	MPP100545	SS-T5A	92713
10-6275	545	MPP100545Z	SS-T5A	92713
10-6280	555	MPP100555	SS-T6	92152
10-6285	555	MPP100555Z	SS-T6	92152
10-6290	558	MPP100558	SS-T7	92318
10-6295	558	MPP100558Z	SS-T7	92318
10-6300	215	MPP100215	3	90408
10-6305	215	MPP100215Z	3	90408
10-6310	225	MPP100225	6	90610
10-6315	225	MPP100225Z	6	90610
10-6320	230	MPP100250	13	90515
10-6325	230	MPP100250Z	13	90515
10-6330	245	MPP100245	7A	90713
10-6335	245	MPP100245Z	7A	90713
10-6340	255	MPP100255	16	91315
10-6345	255	MPP100255Z	16	91315
10-6350	410S	MPP100410S	R40	—
10-6355	415S	MPP100415S	R41	—
10-6360	420S	MPP100420S	R42	—
10-6365	425S	MPP100425S	R44	—



# BIOHAZARD WASTE BAGS

Biohazard Waste Collection Bags are specifically designed to collect, store and transport hazardous waste within a healthcare facility. Made of LLPDE to resist puncture, leaking and abrasion failures, bags meet federal and state regulations for tear and impact resistance. International Biohazard symbol and Bilingual Labeling for easy identification of hazardous waste helps reduce the risk of cross-contamination during infectious waste handling. Offered in various sizes and individually folded for easy dispensing.



Red liners printed with:

- Danger-Peligro
- Infectious Waste
- Basura Infecciosa
- Biohazard

PRODUCT NUMBER	FEATURES	SIZE	GAUGE	CAPACITY	QUANTITY
10-6175	LLPDE, Red w / Black Print	24" x 23"	1.3 MIL	8-10 Gallon	500 Flat / Cs
10-6180	LLPDE, Red w / Black Print	24" x 32"	1.3 MIL	12-16 Gallon	500 Flat / Cs
10-6185	LLPDE, Red w / Black Print	30" x 43"	1.3 MIL	20 Gallon	200 Flat / Cs
10-6190	LLPDE, Red w / Black Print	30" x 43"	3.0 MIL	20 Gallon	100 Flat / Cs
10-6195	LLPDE, Red w / Black Print	33" x 39"	1.3 MIL	33 Gallon	150 Flat / Cs
10-6200	LLPDE, Red w / Black Print	40" x 46"	1.3 MIL	40-45 Gallon	200 Flat / Cs
10-6205	LLPDE, Red w / Black Print	37" x 50"	1.3 MIL	44 Gallon	150 Flat / Cs
10-6210	LLPDE, Red w / Black Print	37" x 50"	3.0 MIL	44 Gallon	75 Flat / Cs

Meets DOT, 49 CFR 173.197 requirements for tear resistance of 480 grams (ASTM D1922-94a) and impact resistance of 165 grams (ASTM D1709-97) as required for transportation of regulated medical waste.

## CROSS REFERENCE

SKLAR	HERITAGE	MEDICAL ACTION	BERRY PLASTICS COVIDIEN	PITT	MEDLINE
10-6175	A4823PR	F116	LH2423HPR	VG2400R	NON122424
10-6180	A4832PR	45-54	LHF2432HPR	V3200R	NON122433
10-6185	A6043PR	F117	LHF3043HPR	VG4300R	NON123143
10-6190	A6043ZR	—	—	VG4303R	—
10-6195	A6639PR	F118	LHF3339HPR	VG3901R	NON203143
10-6200	A8046PR	F135	LHF4046HPR	VG4600R	NON123848
10-6205	A7450PR	F130	LHF3750HPR	VG5001R	—
10-6210	A7450ZR	F169	LHF4046X5PR	VG5000R	NON304046

# LINEN BAGS

Choose from printed and non-printed, color-coded Linen Collection Bags to handle your soiled or infectious linen and laundry for collection, storage and transport. LLPDE bags offer superior leak protection and are available in a wide range of sizes with International Biohazard symbol and Bilingual Labeling for easy identification of infectious waste reducing the risk of cross-contamination during waste handling.

PRODUCT NUMBER	FEATURES	SIZE	GAUGE	CAPACITY	QUANTITY
10-6145	Non-Printed, White	38" x 45"	1.3 MIL	44 Gallon	250 Flat / Cs
10-6150	Non-Printed, Blue	38" x 45"	1.3 MIL	44 Gallon	250 Flat / Cs
10-6155	Printed: Soiled Linen, Blue	30" x 43"	1.3 MIL	20 Gallon	200 Flat / Cs
10-6160	Printed: Soiled Linen, Blue	37" x 50"	1.3 MIL	44 Gallon	150 Flat / Cs
10-6165	Printed: Infectious Linen, Yellow	30" x 43"	1.3 MIL	20 Gallon	200 Flat / Cs
10-6170	Printed: Infectious Linen, Yellow	37" x 50"	1.3 MIL	44 Gallon	150 Flat / Cs

Yellow Linen Bags meet DOT, 49 CFR 173.197 requirements for tear resistance of 480 grams (ASTM D1922-94a) and impact resistance of 165 grams (ASTM D1709-97) as required for transportation of regulated medical waste.

## CROSS REFERENCE

SKLAR	HERITAGE	MEDICAL ACTION	BERRY PLASTICS	COVIDIEN	PITT	MEDLINE
10-6145	—	305	—	—	—	—
10-6150	—	306	—	—	VH4614L	—
10-6155	A6043PX	285	LHF3043HPBL	—	VG4301L	NON022938
10-6160	A7450PX	51-40	—	—	VG5000L	NON022246
10-6165	A6043PY	47-57	LHF3043HPY	—	VG4300Y	NON024042
10-6170	A7450PY	51-45	HH384616PY	—	VG5000Y	NON024046

**DANGER  
INFECTIOUS  
LINEN**



**PELIGRO  
SABANAS  
SUCIAS**

**SOILED LINEN  
SABANAS SUCIAS  
LE LINGE SALE**



# SPECIMEN TRANSPORT BAGS

Offering superior leak resistance, Specimen Transport Bags are suitable for transporting body fluids and tissue and contain an extra pouch for paperwork.

PRODUCT NUMBER	FEATURES	SIZE	GAUGE	QUANTITY
10-6370	Dual Pouch, 3 Wall, Orange / Black	6" x 9"	2 MIL	100 Flat / Bx, 10 Bx / Cs
10-6375	Dual Pouch, 3 Wall, Orange / Black	8" x 10"	2 MIL	100 Flat / Bx, 10 Bx / Cs
10-6380	Dual Pouch, 3 Wall, Red / Stat	6" x 9"	2 MIL	100 Flat / Bx, 10 Bx / Cs

Meets DOT, 49 CFR 173.197 requirements for tear resistance of 480 grams (ASTM D1922-94a) and impact resistance of 165 grams (ASTM D1709-97) as required for transportation of regulated medical waste.

## CROSS REFERENCE

SKLAR	MEDICAL ACTION	MEDLINE
10-6370	4075	ARV9853600
10-6375	4915-ORG	—
10-6380	4001	—



# RECLOSABLE BAGS

Convenient zip-style reclosable bags provide for collection and transport of patient items and come in a wide range of sizes.

PRODUCT NUMBER	FEATURES	SIZE	GAUGE	QUANTITY
10-6385	Single Zipper, Clear	4" x 7"	2 MIL	1000 Flat / Cs
10-6390	Single Zipper, Clear	6" x 6"	2 MIL	1000 Flat / Cs
10-6395	Single Zipper, Clear	6" x 8"	2 MIL	1000 Flat / Cs
10-6400	Single Zipper, Clear	8" x 8"	2 MIL	1000 Flat / Cs
10-6405	Single Zipper, Clear	8" x 10"	2 MIL	1000 Flat / Cs
10-6410	Single Zipper, Clear	10" x 12"	2 MIL	1000 Flat / Cs

## CROSS REFERENCE

SKLAR	HERITAGE	PITT	MEDLINE	MEDICAL ACTION
10-6385	J0706SC	SS-5-C	—	Z2.0406
10-6390	—	—	IVPCV861Z	Z2.0606
10-6395	—	—	—	Z2.0608
10-6400	J0708WC	SS-8-C	—	Z2.0808
10-6405	—	—	—	Z2.0810
10-6410	J1011WC	SS-11-C	—	Z2.1012


# POST MORTEM BAGS

Made of strong, high grade plastic with bound edges and heavy duty, rust-proof nylon zippers. Our designs conform to OSHA's universal precautions regarding containment of body fluids and protection from blood-borne pathogens.

PRODUCT NUMBER	FEATURES	SIZE	GAUGE	QUANTITY
10-6215	Adult w/Tags, HD Straight, Zipper, Latex Free, White	36" x 90"	8 MIL	10 Flat / Cs
10-6220	Pediatric w/Tags, HD Straight, Zipper, Latex Free, White	30" x 48"	8 MIL	30 Flat / Cs
10-6225	Infant w/Tags, HD Straight, Zipper, Latex Free, White	18" x 30"	8 MIL	30 Flat / Cs

## CROSS REFERENCE

SKLAR	MEDLINE	BUSSE	PRECISION DYNAMICS	EM ADAMS	BANTA
10-6215	NON70540W	901	55-7000	1601	950359
10-6220	NON70550W	905	55-7100	1606	—
10-6225	NON70560W	906	—	—	—





# FLUOROSCOPE COVERS

PRODUCT NUMBER	FEATURES	SIZE	QUANTITY
96-5088	Non-Sterile	18" x 20"	100 / Cs
96-5089	Non-Sterile	18" x 20"	50 / Cs
96-5079	Non-Sterile	18" x 24"	100 / Cs
96-5080	Non-Sterile	18" x 24"	50 / Cs
96-5091	Non-Sterile	30" x 36"	100 / Cs
96-5092	Non-Sterile	30" x 36"	50 / Cs
96-1086	Non-Sterile	35" x 43"	100 / Cs
96-1087	Non-Sterile	35" x 43"	50 / Cs
96-5094	Non-Sterile	36" x 30"	100 / Cs
96-5095	Non-Sterile	36" x 30"	50 / Cs
96-5082	Non-Sterile	36" x 44"	100 / Cs
96-5083	Non-Sterile	36" x 44"	50 / Cs
96-5085	Non-Sterile	44" x 60"	100 / Cs
96-5086	Non-Sterile	44" x 60"	50 / Cs
96-5560	Non-Sterile	52" x 52"	100 / Cs
96-5562	Non-Sterile	52" x 52"	50 / Cs
96-5561	Non-Sterile	84" x 110"	100 / Cs
96-5563	Non-Sterile	84" x 110"	50 / Cs

## CROSS REFERENCE

SKLAR	ADVANCE	CARDINAL	CFI	MICROTEK
96-5088	—	—	—	—
96-5089	—	—	—	—
96-5079	—	—	—	—
96-5080	—	—	—	—
96-5091	01-3036	6303RT	CFI-667 NS	3036RT
96-5092	01-3036	6303RT	CFI-667 NS	3036RT
96-1086	—	—	CFI-657 NS	3636RT
96-1087	—	—	CFI-657 NS	3636RT
96-5094	01-3630	—	CFI-655 NS	3630RT
96-5095	01-3630	—	CFI-655 NS	3630RT
96-5082	01-3648	—	—	3648RT
96-5083	01-3648	—	—	3648RT
96-5085	—	—	—	4060RT
96-5086	—	—	—	4060RT
96-5560	01-5050	—	—	5050RT
96-5562	01-5050	—	—	5050RT
96-5561	—	—	—	—
96-5563	—	—	—	—

# SETUP COVERS

Setup covers are useful for isolating tables, trays and other surfaces where adhesion is not needed.

PRODUCT NUMBER	FEATURES	SIZE	QUANTITY
96-5065	Non-Sterile	44" x 70"	100 / Cs
96-5066	Non-Sterile	44" x 70"	50 / Cs
96-5068	Non-Sterile	54" x 80"	100 / Cs
96-5069	Non-Sterile	54" x 80"	50 / Cs
96-5071	Non-Sterile	54" x 90"	100 / Cs
96-5072	Non-Sterile	54" x 90"	50 / Cs

## CROSS REFERENCE

SKLAR	ADVANCE	CFI
96-5065	5470CNS	CFI-702 NS
96-5066	5470CNS	CFI-702 NS
96-5068	—	—
96-5069	—	—
96-5071	5490CNS	—
96-5072	5490CNS	—



# DOME COVERS

PRODUCT NUMBER	FEATURES	SIZE	QUANTITY
96-5050	Non-Sterile	14" Deep	100 / Cs
96-5051	Non-Sterile	14" Deep	50 / Cs
96-1101	Non-Sterile	18" Deep	100 / Cs
96-1103	Non-Sterile	18" Deep	50 / Cs
96-5053	Non-Sterile	22" Deep	100 / Cs
96-5054	Non-Sterile	22" Deep	50 / Cs
96-5056	Non-Sterile	27" Deep	100 / Cs
96-5057	Non-Sterile	27" Deep	50 / Cs
96-1109	Non-Sterile	30" Deep	100 / Cs
96-1104	Non-Sterile	30" Deep	50 / Cs

## CROSS REFERENCE

SKLAR	3M	ADVANCE	CFI	MICROTEK
96-5050	—	03-KP14NS	CFI-763NS	0030S
96-5051	—	03-KP14NS	CFI-763NS	0030S
96-1101	1096NS	03-KP18NS	CFI-911NS	0036S
96-1103	1096NS	03-KP18NS	CFI-911NS	0036S
96-5053	1097NS	03-KP22NS	CFI-821NS	0040S
96-5054	1097NS	03-KP22NS	CFI-821NS	0040S
96-5056	—	03-KP26NS	CFI-756NS	0052S
96-5057	—	03-KP26NS	CFI-756NS	0052S
96-1109	—	03-KP30NS	CFI-785NS	0060S
96-1104	—	03-KP30NS	CFI-785NS	0060S

# MAYO STAND COVERS

PRODUCT NUMBER	FEATURES	SIZE	QUANTITY
96-5566	Non-Sterile	20" x 54"	100 / Cs
96-5570	Non-Sterile	20" x 54"	50 / Cs
96-5135	Non-Sterile	23" x 54"	100 / Cs
96-5136	Non-Sterile	23" x 54"	50 / Cs
96-5567	Non-Sterile	30" x 57"	100 / Cs
96-5571	Non-Sterile	30" x 57"	50 / Cs

## CROSS REFERENCE

SKLAR	3M	CARDINAL	CFI	KIMBERLY CLARK	MEDLINE
96-5566	—	—	—	—	—
96-5570	—	—	—	—	—
96-5135	9036NS	A8337N	CFI-708 NS	79601	DYNJP2500
96-5136	9036NS	A8337N	CFI-708 NS	79601	DYNJP2500
96-5567	—	A8339N	—	79602	DYNJP2510
96-5571	—	A8339N	—	79602	DYNJP2510



# BAND BAGS WITH TAPE

PRODUCT NUMBER	FEATURES	SIZE	QUANTITY
96-5000	Non-Sterile	8" x 8"	100 / Cs
96-1059	Non-Sterile	10" x 10"	100 / Cs
96-5582	Non-Sterile	12" x 12"	100 / Cs
96-5001	Non-Sterile	12" x 24"	100 / Cs
96-5515	Non-Sterile	15" x 15"	100 / Cs
96-5516	Non-Sterile	20" x 18"	100 / Cs
96-5583	Non-Sterile	20" x 20"	100 / Cs
96-1061	Non-Sterile	20" x 28"	100 / Cs
96-5002	Non-Sterile	30" x 15"	100 / Cs
96-1062	Non-Sterile	30" x 20"	100 / Cs
96-5517	Non-Sterile	30" x 30"	100 / Cs
96-5584	Non-Sterile	30" x 36"	100 / Cs
96-5518	Non-Sterile	36" x 15"	100 / Cs
96-5010	Non-Sterile	36" x 20"	100 / Cs
96-5580	Non-Sterile	36" x 24"	100 / Cs
96-1068	Non-Sterile	36" x 28"	100 / Cs
96-5581	Non-Sterile	36" x 30"	100 / Cs
96-5004	Non-Sterile	36" x 36"	100 / Cs
96-5585	Non-Sterile	36" x 40"	100 / Cs
96-5579	Non-Sterile	36" x 48"	100 / Cs
96-5586	Non-Sterile	36" x 50"	100 / Cs
96-1052	Non-Sterile	36" x 54"	100 / Cs
96-5529	Non-Sterile	40" x 30"	100 / Cs
96-5530	Non-Sterile	40" x 40"	100 / Cs
96-5018	Non-Sterile	40" x 60"	100 / Cs
96-5007	Non-Sterile	42" x 36"	100 / Cs
96-5531	Non-Sterile	44" x 36"	100 / Cs
96-1047	Non-Sterile	44" x 40"	100 / Cs
96-5008	Non-Sterile	50" x 36"	100 / Cs
96-5012	Non-Sterile	50" x 40"	100 / Cs
96-5009	Non-Sterile	50" x 50"	100 / Cs
96-5604	Non-Sterile	54" x 90"	100 / Cs
96-5600	Non-Sterile	60" x 30"	100 / Cs

## CROSS REFERENCE

SKLAR	ADVANCE	CARDINAL	CFI	MICROTEK
96-5000	—	—	—	—
96-1059	01-1010	—	—	1010RT
96-5582	01-1212	—	—	1212RT
96-5001	—	—	—	1224RT
96-5515	—	—	—	1515RT
96-5516	—	—	CFI-664 NS	2018RT
96-5583	01-2020	—	CFI-651 NS	2020RT
96-1061	01-2028	—	—	2028RT
96-5002	—	—	—	3015RT
96-1062	01-3020	—	—	3020RT
96-5517	01-3030	30309RT	CFI-661 NS	3030RT
96-5584	01-3036	6303RT	CFI-667 NS	3036RT
96-5518	—	3620RT	CFI-833 NS	3625RT
96-5010	01-3620	—	CFI-654 NS	3620RT
96-5580	—	—	CFI-932 NS	3624RT
96-1068	01-3628	—	CFI-835 NS	3628RT
96-5581	01-3630	—	CFI-655 NS	3630RT
96-5004	01-3636	—	CFI-657 NS	3636RT
96-5585	01-3640	—	CFI-879 NS	3640RT
96-5579	01-3648	—	—	3648RT
96-5586	—	—	—	3650RT
96-1052	01-3654	—	—	3654RT
96-5529	01-4030	—	CFI-761 NS	4030RT
96-5530	01-4040	—	CFI-827 NS	4040RT
96-5018	01-4060	—	—	4060RT
96-5007	—	—	—	4236RT
96-5531	—	—	CFI-777T NS	4436RT
96-1047	—	—	—	4440RT
96-5008	—	—	—	5036RT
96-5012	01-4440	—	—	5040RT
96-5009	01-5050	—	—	5050RT
96-5604	—	—	—	5490RT
96-5600	—	—	—	6030RT

SKLAR brand banded bags are a multipurpose covering with an elastic opening. Quick and easy to fit over a variety of OR equipment.



# WASTE CAN LINERS

Designed for non-medical waste disposal, Waste Can Liners are used throughout health care facilities to collect and contain refuse and are available in a wide variety of sizes and film gauges to meet your specific application. Clear liners allow for quick visualization of contents before disposal while colored options hide unsightly waste products.

Choose the right plastic types:

- **LLDPE (Linear Low Density Polyethylene)** provides maximum puncture and tear resistance under tough transport conditions and is best suited for handling objects with sharper edges.
- **HDPE (High Density Polyethylene)** provides high resistance to puncturing, durability and greater load capacity for paper or non-sharp objects under moderate transport conditions.

All bags are designed with Star Seals to eliminate gaps, maximize carrying capacity and minimize leaks and conveniently offered in flat packs or coreless rolls for easy dispensing.

PRODUCT NUMBER	FEATURES	SIZE	GAUGE	CAPACITY	QUANTITY
10—6000	HDPE, Natural	24" x 33"	8 MIC	12-16 Gallon	50/Roll, 20 Rolls/Cs
10—6005	HDPE, Natural	30" x 37"	10 MIC	20-30 Gallon	25/Roll, 20 Rolls/Cs
10—6010	HDPE, Natural	30" x 37"	10 MIC	20-30 Gallon	500 Flat / Cs
10—6015	HDPE, Natural	30" x 37"	16 MIC	20-30 Gallon	25/Roll, 20 Rolls/Cs
10—6020	HDPE, Natural	30" x 37"	16 MIC	20-30 Gallon	500 Flat / Cs
10—6025	HDPE, Natural	33" x 40"	12 MIC	33 Gallon	25/Roll, 10 Rolls/Cs
10—6030	HDPE, Natural	33" x 40"	12 MIC	33 Gallon	250 Flat / Cs
10—6035	HDPE, Natural	33" x 40"	16 MIC	33 Gallon	25/Roll, 20 Rolls/Cs
10—6040	HDPE, Natural	33" x 40"	16 MIC	33 Gallon	500 Flat / Cs
10—6045	HDPE, Natural	40" x 48"	13 MIC	40-45 Gallon	25/Roll, 10 Rolls/Cs
10—6050	HDPE, Natural	40" x 48"	13 MIC	40-45 Gallon	250 Flat / Cs
10—6055	HDPE, Natural	40" x 48"	16 MIC	40-45 Gallon	25/Roll, 10 Rolls/Cs
10—6060	HDPE, Natural	40" x 48"	17 MIC	40-45 Gallon	25/Roll, 10 Rolls/Cs
10—6065	HDPE, Natural	40" x 48"	17 MIC	40-45 Gallon	250 Flat / Cs
10—6070	HDPE, Natural	43" x 48"	16 MIC	60 Gallon	25/Roll, 8 Rolls/Cs
10—6075	LLDPE, Clear	24" x 32"	0.6 MIL	12-16 Gallon	500 Flat / Cs
10—6080	LLDPE, Clear	30" x 36"	0.8 MIL	20-30 Gallon	250 Flat / Cs
10—6085	LLDPE, Clear	33" x 39"	0.6 MIL	33 Gallon	150 Flat / Cs
10—6090	LLDPE, Clear	40" x 46"	0.8 MIL	40-45 Gallon	125 Flat / Cs
10—6095	LLDPE, Clear	38" x 45"	1.0 MIL	44 Gallon	125 Flat / Cs
10—6100	LLDPE, Black	43" x 48"	1.0 MIL	56 Gallon	100 Flat / Cs
10—6105	LLDPE, Super Heavy, Clear	33" x 39"	1.1 MIL	33 Gallon	250 Flat / Cs
10—6110	LLDPE, Super Heavy, Gray	33" x 39"	1.1 MIL	33 Gallon	250 Flat / Cs
10—6115	LLDPE, Super Heavy, Clear	40" x 46"	1.1 MIL	40-45 Gallon	125 Flat / Cs
10—6120	LLDPE, Super Heavy, Gray	40" x 46"	1.1 MIL	40-45 Gallon	125 Flat / Cs
10—6125	LLDPE, Super Heavy, Clear	43" x 47"	1.1 MIL	56 Gallon	100 Flat / Cs
10—6130	LLDPE, Super Heavy, Gray	43" x 47"	1.1 MIL	56 Gallon	100 Flat / Cs
10—6135	LLDPE, Super Heavy, Clear	38" x 55"	1.1 MIL	60 Gallon	100 Flat / Cs
10—6140	LLDPE, Super Heavy, Gray	38" x 55"	1.1 MIL	60 Gallon	100 Flat / Cs

## CROSS REFERENCE

SKLAR	HERITAGE	MEDICAL ACTION	BERRY PLASTICS COVIDIEN	PITT	MEDLINE
10—6000	Z4833LN R01	RS243308N	HR243308N	VVMR24330MC	NONCHD33MC
10—6005	Z6037MN R01	RS303710N	HR303710N	VMR30371MC	NONCHD3710C
10—6010	Z6037MN	HCL303710	—	VH30371EC	NONCHD3710C
10—6015	Z6037XN R01	RS303716N	HR303716N	VMR30373MC	NONCHD37MC
10—6020	Z6037XN	—	—	VH30373EC	NONCHD37MC
10—6025	Z6640MN R01	RS334011N	HH334013N	VMR33402MC	NONCHD40HC
10—6030	Z6640MN	HCL334011	—	H33402EC	NONCHD40HC
10—6035	Z6640VN R01	RS334016N	HR334016N	VMR33403MC	NONCHD40XC
10—6040	Z6640VN	HCL334016	—	VH33403EC	NONCHD40XC
10—6045	Z8048HN R01	RS404812N	HR404812N	VMR40482MC	NONCHD48HC
10—6050	Z8048HN	HCL404812	—	VH40482EC	NONCHD48HC
10—6055	Z8048VN R01	RS404816N	HR404816N	VMR40483MC	NONCHD48XC
10—6060	Z8048XN R01	RS404817N	—	VMR40483MC	—
10—6065	Z8048XN	—	—	VH40483EC	—
10—6070	Z7660VN R01	RS434816N	HR434816N	VMR38603MC	NONCHD47HC
10—6075	D4832RC	X2418	LSF2325LC	VP3310C	NONO222207
10—6080	H6036MC	4637	LHF3036XR	VP3715C	NONO22220C
10—6085	H6639TC	R4646, 4644	LHF3340HC	VP4015C	NONO222228
10—6090	H8046HC	4660	LSF4046XC	VP4815C	NONO22232C
10—6095	H7444HC	4659	—	—	—
10—6100	H8647HB	X2414	—	VP4940C	—
10—6105	H6639SC	1286	LPF3340X3G	VP4030C	—
10—6110	H6639SG	4648	LPF3340X3C	VP4830C	—
10—6115	H8046SC	23-50C	LPF4046X3G	VRX483XG	—
10—6120	H8046SG	3916	LPF4046X3C	VP4830C	—
10—6125	H8647SC	4666	LPF4347X3G	VP4930C	—
10—6130	H8647SG	4665	LPF4347X3C	VP4730C	—
10—6135	H7658SC	23-65C	LHF3858HPR	VP6530C	—
10—6140	H7658SG	4661	—	VB76530K	—





READY-TO-USE, PRE-MOISTENED TOWELETTES

# DISINFECTANT™ SURFACE WIPES



## DISINFECTANT™ SURFACE WIPES

**2-Minute Kill Time!**

**SAME EFFECTIVE FORMULA AS  
SKLAR DISINFECTANT!**

- Ready to use, non-toxic formula
- Specially formulated for the control of cross contamination
- Effective against a broad range of pathogens including H1N1, HIV-1, HPV, VRE, Norovirus, Hep B & C as well as TB, MRSA, A. Baumannii, and C. Diff (veg.).
- Safe for all areas of the facility including the NICU, operating rooms, isolation rooms, patient care areas and laboratories
- Satisfies OSHA's Bloodborne Pathogen Standards
- Destroys odor-causing bacteria
- 7" x 10" non-woven, tear resistant towelette
- Clean Scent
- EPA registered

### PRODUCT NUMBER

**10-1616**

FEATURES: 100 TOWELETTES PER CANISTER

QUANTITY: 6 / CASE

### PRODUCT NUMBER

**10-1617**

FEATURES: 100 TOWELETTES PER CANISTER

QUANTITY : 12 / CASE

CONVENIENT  
DISPENSING  
CANISTER

**KILLS MRSA AND H1N1 VIRUS  
2-MINUTE KILL TIME!**

Please direct your questions or comments about these products to:



889 South Matlack Street  
West Chester, PA 19382

Telephone: 800.221.2166 / 610.430.3200

Corporate Fax: 610.429.0500

Sales Fax: 610.696.9007

Email: [international@sklarcorp.com](mailto:international@sklarcorp.com)

[surgi@sklarcorp.com](mailto:surgi@sklarcorp.com)

[www.sklarcorp.com](http://www.sklarcorp.com)

**Distributed By:**